

all-in-ONE™

Surgical Instrument Cleaner AND Conditioners

Highly concentrated for fast and effective cleaning.

Image: all-in-ONE™ Surgical Instrument Cleaner AND Conditioners replacing multiple products, lowering costs, and improving the performance of Surgical Instruments. The all-in-ONE™ Cleaners AND Conditioners deliver a unique formulation of detergents and multi-tiered "high-level" enzymes (Protease, Amylase, Lipase, and Carbohydrase) within a Neutral pH surfactant chemical complex. The all-in-ONE™ rapidly breaks down and removes all forms of proteinaceous bioburden, cleans and conditions the surface, enhances the "passive layer" protection of surgical stainless steel, lubricates moving parts, and then rinses clean. **ONE Gallon** of the all-in-ONE™ highly concentrated Surgical Instrument Cleaner AND Conditioners will deliver **512 gallons** (1938 liters) of Surgical Instrument Cleaning AND Conditioning power. The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and **STAINS**.

The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and **STAINS**.

The **easy FOAM-it™** lays down a thick blanket to hydrate and **remove** proteinaceous debris, hard water deposits, and encrusted **STAINS**, for fast acting enzyme detergent cleaning action. When water is added the 'foam' dissipates to provide a **clear** 'free-rinsing' cleaning AND conditioning solution.

all-in-ONE™ Surgical Instrument Cleaner AND Conditioners

ONE product
for cleaning and conditioning Surgical Instruments.
ONE application
to thoroughly clean and condition Surgical Instruments.
Apply it **ONCE...**
It's all you really need for cleaning and conditioning
Surgical Instruments and Scopes.

Replace Cleaning AND Conditioning products for:
removing all forms of bioburden,
pre-cleaning and soaking to remove **STAINS**,
manual cleaning, Ultrasonic cleaning,
and Washer Decontaminator Disinfectant reprocessing.

The all-in-ONE™ Surgical Instrument Cleaners deliver a unique formulation of detergents and multi-tiered "high-level" enzymes (Protease, Amylase, Lipase, Carbohydrase) within a Neutral pH surfactant chemical complex. The all-in-ONE™ cleaners will: break down bioburden and clean the surface, remove **STAINS**, enhance the "passive layer" protection of surgical stainless steel, and lubricate moving parts. The highly concentrated all-in-ONE™ cleaners are aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment.

easy FOAM-it™ Surgical Instrument Cleaning AND Conditioning 12 applicator bottles per case

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 4 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 5 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 15 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 30 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 55 gallons

Your cleaner surgical instruments will look newer and perform longer.

yourCEBA

Customer Service: 1-916-335-6113 Technical Support: 509-747-5027

Website: www.yourCEBA.com, Email: Support@yourCEBA.com, Email: Info@yourCEBA.com



easy FOAM-it™
Surgical Instrument
Cleaners AND Conditioners



all-in-ONE™
Surgical Instrument
Cleaners AND Conditioners